PD Percept Pipeline

Samantha N. Pitts – Edited (November 2023)

Shahed Lab

Neurology, Mount Sinai West/ Icahn School of Medicine


Samantha.pitts@mssm.edu

Version1 (Published): https://github.com/C-ACT/Percept-LFP-Toolbox/

Version2: Percept Extract Code

Scripts: Working Upload	Version of Matlab Used: 2021a - 2022b
	Prompts a single BrainSense Medtronic JSON file from up to 60 days prior to office visit. Extracts, saves and plots
	indefinite streaming and timeline data and stimulation group settings in folder of choice as "LFPTrendlogs.mat" structure and figure.
loadJSONs	As adapted by Yohann Thenaisie 02.09.2020 - Lausanne University Hospital (CHUV)
continue Course	Synthesises a cumulative dataset (as a .mat file) of stimulation settings from "LFPTrendlogs.mat <activegroup.": a="" and="" create="" extracted="" from="" jsons="" lfptrendlogs="" loads="" of="" patient,="" prompts="" queue="" select="" single="" td="" the<="" to="" user=""></activegroup.":>
combine_Groups	Active Group Variable into singular table, and saves as "PD_AG.mat." Repetitive data is overwritten.
structureData	Generates a cumulative dataset (as a .mat file) of stimulation amplitudes and raw LFP band power from "LFPTrendlogs.mat <lfptrendlogs" "dates",="" "hh:mm:ss"="" "structuredata.mat"="" (default="" (unknown="").="" ,="" a="" added="" and="" are="" as="" assigned="" bounds,="" but="" c1="09:00:00" c2="21:00:00" c2.="" chronologically,="" constraint="" contrained="" data="" data.="" dataset="" duplicate="" each="" extracted="" extraction="" fill="" format,="" from="" generated="" in="" is="" jsons="" lfptrendlogs="" lower="" missing="" nan="" of="" ommission="" patient.="" prompts="" ranges="" reflects="" respectively.="" result="" row="" s="" saved="" select="" set="" single="" sorted="" sourced="" specified="" struct="" t="" table="" tables="" td="" the="" time="" times="" to="" upper="" user="" value)<="" variables="" via="" with=""></lfptrendlogs">
visData	Calculates the daily IQR, daily standard deviation and FS using a non-overlapping window of k, indicated range of days (default is k=5). Prompts user to select "structureData.mat," or cumulative JSON data according to timetable T and S. @igr matlab function is applied for each unique date of data, across each hemisphere LFP band power row within the T timetable, then is @std. Fluctuation score is calculated after filling each null value igr row (date/day) with the previous known value. k consecutively summed values of IQR, are reported as day 1 out of k FS. Process is repeated until the last group of k days has its value recorded. These statistics are saved to folder. Standard deviation and fluctuation scores are plotted
saveCSV	Timetables are saved as csv
PD_Timeline_Extended (pd_01)	Prompts for table that includes the date of first programming. Reorganizes and rereferences data according to
Session Category	provided first programming date and subsequent visits or "sessions". Puts all measurements sequentially in a
StackTablesVariables	single column for JMP compatibility/ analysis.
Resampling and Checking Groups	
oadJSONs2/ExtractGroups	Re-extracts stimulation settings to accommodate found inconsistencies (missing samples, sample hemisphere
combineAllGroups	allocation) and inclusions of more precise setting changes (pulse width, rate, band of interest, contact)

Script: LoadJSONs



91649 struct

LFPTrendLogs 1x1



Contents of the Active Group Variable contains a table of the Stimulation settings and "Groups" Streamed by user.



Contents of the LFPTrendLogs Variable contains a struct of stimulation amplitude data and LFP Band Power averaged sampled per 10 minutes (REF). Additional contents unused include user-recorded events.


Script: Combine_Groups

	DateTime	OldGroupId	NewGroupId	NewProgramSettings
1	03-Nov-2022 20:45:50	'GROUP_A'	'GROUP_A'	1×1 struct
2	19-Oct-2022 19:23:23	'GROUP_A'	'GROUP_A'	1×1 struct
AG:	2×4 table			
	DateTime	OldGroupId	NewGroupId	NewProgramSettings
1	19-Oct-2022 19:23:23	'GROUP_A'	'GROUP_A'	1×1 struct
2	19-Oct-2022 19:23:23	'GROUP A'	'GROUP A'	1×1 struct
	2×4 table	GROOF_A	GROUP_A	r~r struct
		OldGroupId	NewGroupId	
	2×4 table			
AG :	= 2×4 table DateTime	OldGroupId	NewGroupId	NewProgramSettings
AG :	2×4 table DateTime 19-Oct-2022 19:23:23	OldGroupId 'GROUP_A'	NewGroupId 'GROUP_A'	NewProgramSettings
AG :	DateTime 19-Oct-2022 19:23:23 19-Oct-2022 19:23:23	OldGroupId 'GROUP_A'	NewGroupId 'GROUP_A'	NewProgramSetting: 1×1 struct 1×1 struct
AG :	= 2×4 table DateTime 19-Oct-2022 19:23:23 19-Oct-2022 19:23:23 2×4 table	OldGroupId 'GROUP_A' 'GROUP_A'	NewGroupId 'GROUP_A' 'GROUP_A'	NewProgramSettings

PD_AG = 2×3 timetable										
		DateTime	OldGroupId	NewGroupId	NewProgramSettings					
1	1	19-Oct-2022 19:23:23	'GROUP_A'	'GROUP_A'	1×1 struct					
2	2	13-Jan-2023 15:31:49	'GROUP_A'	'GROUP_A'	1×1 struct					

Script: structureData

Each Table contains an array of data values, specified by time corresponding to each date column.




	1	1 2		4	5	6
Date	Raw Date Time (HH	Raw LFP Band Powe	Duration (HH:mm:s	LFP Band Power - R/CH1	LFP Band Power	Date Time (HH:mm:ss
05-May-2021	1x53 datetime	53x2 double	1x35 duration	35x1 double	0	1x35 datetime
06-May-2021	1x144 datetime	144x2 double	1x72 duration	72x1 double	0	1x72 datetime
07-May-2021	1x144 datetime	144x2 double	1x72 duration	72x1 double	0	1x72 datetime
08-May-2021	1x144 datetime	144x2 double	1x72 duration	72x1 double	0	1x72 datetime
09-May-2021	1x141 datetime	141x2 double	1x72 duration	72x1 double	0	1x72 datetime
10-May-2021	1x112 datetime	112x2 double	1x59 duration	59x1 double	0	1x59 datetime
11-May-2021	144x1 datetime	144x2 double	72x1 duration	72x1 double	0	72x1 datetime
12-May-2021	144x1 datetime	144x2 double	72x1 duration	72x1 double	0	72x1 datetime
12 May 2021	1 4 4 1	144.2 4	721 d	721	n	72.1 detetions

Table T. Note Fluctuation score provided is not yet log'd.


	1	2	: 3	4	5	6	7
Date	Duration (HH:mm:ss)	Stim Amp - R/CH1 (mA)	Stim Amp - L	Channel Names	Band of Intere	Pulsew	Date Time (HH:n
05-May-2021	1x107 duration	107x1 double		"ONE_AND_THREE"	15.6300	60	1x107 datetime
06-May-2021	1x144 duration	144x1 double		"ONE_AND_THREE"	15.6300	60	1x144 datetime
07-May-2021	1x216 duration	216x1 double	0	"ONE_AND_THREE"	15.6300	60	1x216 datetime
08-May-2021	1x216 duration	216x1 double	0	"ONE_AND_THREE"	15.6300	60	1x216 datetime
09-May-2021	1x216 duration	216x1 double		"ONE_AND_THREE"	15.6300	60	1x216 datetime
10-May-2021	1x203 duration	203x1 double	0	"ONE_AND_THREE"	15.6300	60	1x203 datetime
11-May-2021	72x1 duration	72x1 double	0	"ONE_AND_THREE"	17.5800	60	72x1 datetime
12-May-2021	72x1 duration	72x1 double	0	"ONE_AND_THREE"	17.5800	60	72x1 datetime
13_May_2021	72v1 duration	72v1 double	п	"ONE AND THREE"	17 5800	60	72v1 datetime
		7	able S				

Script: visData




Calculated outputs, (standard deviation and fluctuation score). *Note fluctuation score must be log(x) d in data tables/raw tables, but depiction is a log-x-axis scale.

Old version of this script (Github) provides the following image sequences:



Additional Examples



Script: PD_Timeline_Extended (pd_01)

PD_Timeline_Extended								
Days Since First Programming	Data	Group	Patient	Data_Type	Session Number			
-22	NaN	GROUP_B		IQR (uVp) - Left	1			
-22	NaN	GROUP_B		IQR (uVp) - Right	1			
-21.6206481481481	0	GROUP_B		Stimulation Amplitude (mA) - Left	1			
-21.6206481481481	0	GROUP_B		Stimulation Amplitude (mA) - Right	1			
-21.6137037037037	0	GROUP_B		Stimulation Amplitude (mA) - Left	1			
-21.6137037037037	0	GROUP_B		Stimulation Amplitude (mA) - Right	1			
46.7057523148148	0	GROUP_B		Stimulation Amplitude (mA) - Left	2			
46.7057523148148	0	GROUP_B		Stimulation Amplitude (mA) - Right	2			
46.7057523148148	0	GROUP_B		Stimulation Amplitude (mA) - Right	2			
46.7093634259259	1.7	GROUP_B		Stimulation Amplitude (mA) - Left	2			
46.7093634259259	1.9	GROUP_B		Stimulation Amplitude (mA) - Right	2			
46.7093634259259	83	GROUP_B		LFP Band Power - Left	2			
46.7093634259259	286	GROUP_B		LFP Band Power - Right	2			
46.7126967592593	0	GROUP_B		Stimulation Amplitude (mA) - Left	2			

Limitations/Improvements from first Version (Github)

- In PD_AG Generation the following needs to be automated
- Rename groups if the, BOI, and channel is different only
- A Question if we group the left and right as different classified "groups" even if one hemisphere remains the same? Maybe this explains variation later so for now separate left and right hemisphere groups
- Case of Different Group, Same Settings
- Normalize LFP to reduce outliers before FS calculated

Data Statistics Steps in JMP (V.2022-23)

Recording and re-referencing according to program dates (JSON file)

Each JSON file corresponds to an office visit, which is notated below. Data used to separate data acquired between office visits by assigning data the *session number* category variable.

RedCap									
Subject	Record ID	1st programming	2nd	3rd	4th	5th	6th		
1	H	11/25/22	12/9/22	1/11/23	2/8/23	3/8/23			
2		6/8/22	7/13/22	8/10/22	10/12/22	12/14/22			
3		7/6/22	8/3/22	10/12/22	11/9/22	2/8/23			
4		5/11/22	5/11/22	6/15/22	7/1/22	8/3/22	8/24/22	9/21/22	
5	FEET TO SERVICE STATE OF THE S	6/8/22	7/8/22	9/7/22	10/5/22	12/7/22	3/10/23		


Pasting in Table

Outcome of last Matlab script

PD_Timeline_Extended

Days Since First Programming	Data	Group	Patient	Data_Type	Session Number
-19	10.5	GROUP_A		IQR (uVp) - Right	1
-19	90.5	GROUP_A		Fluctuation Score (uVp) - Right	1
-19	1216.5	GROUP_A		IQR (uVp) - Left	1
-19	6647	GROUP_A		Fluctuation Score (uVp) - Left	1
-18.205787037037	0	GROUP_A		Stimulation Amplitude (mA) - Left	1
-18.205787037037	0	GROUP_A		Stimulation Amplitude (mA) - Right	1
-18.205787037037	89	GROUP_A		LFP Band Power - Right	1
-18.205787037037	2180	GROUP_A		LFP Band Power - Left	1
-18.1988425925926	0	GROUP_A		Stimulation Amplitude (mA) - Left	1
-18.1988425925926	0	GROUP_A		Stimulation Amplitude (mA) - Right	1
-18.1988425925926	86	GROUP_A		LFP Band Power - Right	1
-18.1988425925926	2420	GROUP_A		LFP Band Power - Left	1
-18.1918981481481	0	GROUP_A		Stimulation Amplitude (mA) - Left	1
-18.1918981481481	0	GROUP_A		Stimulation Amplitude (mA) - Right	1
-18.1918981481481	88	GROUP_A		LFP Band Power - Right	1
-18.1918981481481	2098	GROUP_A		LFP Band Power - Left	1
-18.1849537037037	0	GROUP_A		Stimulation Amplitude (mA) - Left	1
-18 1849537037037		GROUP A		Stimulation Amplitude (mA) - Right	1

Paste all patients into single JMP table




*Note: Additional Variable and data type changes. *Session number* was set to numeric ordinal variable. *Days Since first programming* was copied to a new column and set to numeric ordinal; formula applied was as follows to categorize data into days and visualize the distribution of data per day.




- 1. Data Table of all patients, Rows>DataFilter>: greater than -0.0001 and less than 180.0001: Image basics. Will Need to scale images individually/ view at 1-2 levels at a time.
- 2. Manually Filter Table to reflect stimulation settings, including listing all stimulation changes from each json, re-referencing setting start dates, re-grouping settings in categories that reflect unique band power sensing bands and sensing channel/contact, as well as pulse width. Table reordered and referenced to group dataset in JMP by recording setting specifications. EX:


$\overline{}$	Days Elaps ▼	Reassigned Group - L	Reassigned Group - R	GroupNam ▼	Dates	₩	BP - Left ▼	BP - Right	Channel - L	Channel - F PW-L	▼ PW-R
L	-36.1920949	1	1	Α	'19-Oct-2022 19:23:23'		15.63	17.58	ZERO_AND_T	ZERO_AND_TWO	
L	-21.134838	2	2	Α	'03-Nov-2022 20:45:50'		12.7	9.77	ZERO_AND_T	ZERO_AND_TWO	
L	49.64709491	3	3	Α	'13-Jan-2023 15:31:49'		14.65	19.53	ZERO_AND_T	ONE_AND_THREE	
L	-67.157419	1	1	Α	'01-Apr-2022 20:13:19'		14.65	13.67	ZERO_AND_T	ZERO_AND_THREE	
L	0.570162037	1	1	Α	'08-Jun-2022 13:41:02'		14.65	13.67	ZERO_AND_T	ZERO_AND_THREE	
L	0.574421296	2	2	В	'08-Jun-2022 13:47:10'		22.46	62.5	ZERO_AND_T	ZERO_AND_TWO	
Ш	63.58305556	2	3	В	'10-Aug-2022 13:59:36'		22.46	21.48	ZERO_AND_T	ZERO_AND_TWO	
Ш	-8.42806713	1	1	Α	'27-Jun-2022 13:43:35'		21.48	21.48	ZERO_AND_T	ZERO_AND_TWO	
	28.57388889	2	2	Α	'03-Aug-2022 13:46:24'		16.6	23.44	ONE_AND_TH	ONE_AND_THREE	
L	98.58853009	3	3	Α	'12-Oct-2022 14:07:29'		22.46	20.51	ONE_AND_TH	ONE_AND_THREE	
L	-221.095694	4	2	Α	'01-Oct-2021 21:42:12'		83.98	28.32	ZERO_AND_T	ZERO_AND_TWO	
L	-18.2193171	5	4	Α	'22-Apr-2022 18:44:11'		18.55	11.72	ZERO_AND_T	ZERO_AND_THREE	
L	0.598171296	5	4	Α	'11-May-2022 14:21:22'		18.55	11.72	ZERO_AND_T	ZERO_AND_THREE	
L	0.60287037	6	1	В	'11-May-2022 14:28:08'		24.41	24.41	ZERO_AND_T	ONE_AND_THREE	
	35.63908565						(Delete)				

1. Subset: (By Patient, for each patient)






- 1. Split: and Rename by Group Patient Type
- 1. Filter Again Rows>DataFilter>: greater than -0.0001 and less than 180.0001:




1. Exclude data again for statistical testing: Rows>DataFilter>: greater than - 0.0001 and less than 180.0001


2. Interpolating the stimulation amplitudes for missing values; Fill missing Values, interpolate



4. Split Tables into "Left and Right Hemispheres" via Tables>Subset



5. Program Settings Encoded to Settings-R or Settings-L



Filter out beta bands <12.5 >30 Hz

6. Scatterplot Matrix via Analysis>Multivariate

